

Reference Photos for Blood Draw Packing and Shipping to University of California Davis

Supplies: insulated container, bubble wrap, Ziplock bags.
It's important to wrap vial in bubble wrap to prevent damage.
Place the submission form into a separate Ziploc bag and place inside the insulated container.



Size comparison, vials are tiny



FedEx has a variety of small packaging supplies, free of charge.
Boxes and padded envelopes.



Sample of FedEx shipping label to UCD.

From: (530) 752-4892 Stem Lab UC Davis 2108 Tupper Hall One Shields Ave Davis, CA 95616	Origin ID: SMFA		Ship Date: 15AUG18 Act/Wgt: 2.0 LB CAD: 454814FXRS1409
SHIP TO: (530) 752-5058 Amino Acid Laboratory University of California Davis 1020 Vet Med 3B 1089 Veterinary Medicine Drive Davis, CA 95616		BILL SENDER	Delivery Address Bar Code 
			Ref # F3V035150 Invoice # PO # Dept # EEAARNZ
			RELEASE#: 3785346
			THU - 16 AUG 10:30A PRIORITY OVERNIGHT
		TRK# 4420 2401 2925 6281	95616 CA-US SMF
		84 SMFA	
			522C15309F220